Implant Center

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ENGLISH

CONTENTS

1- INTRODUCTION	5	Micro-motor cord and scaler cord maintenance	19
2- WARNINGS	5-7	I-SURGE micro-motor maintenance	19
3- DESCRIPTION		Ultrasonic handpiece maintenance	20
Physical description	7	Tip maintenance	20
Technical description	7-10	12- MONITORING/PREVENTIVE AND CORRECTIVE MAINTENANCE	
4- INSTALLATION/FIRST USE		Monitoring	20
Unpacking the device	10	Preventive and corrective maintenance	21
Recommendations	10	Fuse replacement	21
Installation	10-11	Operation faults	21-22
First use	11	13- ELECTROMAGNETIC COMPATIBILITY	
5- IRRIGATION SOLUTION	11	Electromagnetic emissions	23
6- IRRIGATION LINES	11	Electromagnetic immunity	24
7- ADJUSTMENTS/SETTINGS/MODES/INTERFACE		Electromagnetic immunity / radiofrequency mobile equipment	24-25
Startup	11	Recommended separation distances	25
Settings	11	Cable lengths	26
Irrigation function	11-12	14- DISPOSAL AND RECYCLING	26
I-SURGE function	12-13	15- MANUFACTURER'S LIABILITY	26
PIEZOTOME function	13-14	16- REGULATIONS	27
NEWTRON function	14	17- SYMBOLS	27
Toolbox	14-15	18- CUSTOMER RELATIONS	
8- SAFETY	15	Manufacturer identification	148
9- USING THE DEVICE	16-17	Subsidiaries	148-149
10- SHUTTING DOWN THE DEVICE	17	The list of accessories referred to in this manual is non exhaustive. Consult your local dealer for further details.	
11- ROUTINE MAINTENANCE/STERILIZATION			
Irrigation line maintenance	18		
Contra-angle maintenance	18		
Rotary instrument maintenance	18		
Device maintenance	18-19		

I - INTRODUCTION

Congratulations! You have just taken possession of your IMPLANT CENTER®.

Designed by SATELEC®, the IMPLANT CENTER is a multifunction device intended for:

- dental implantology, with its I-SURGE micro-motor,
- surgical procedures (osteotomy, osteoplasty, periodontal and implant surgery) with its PIEZOTOME surgical ultrasonic handpiece,
- mechanical ultrasonic treatment (prophylaxis, periodontics, endodontics) with its NEWTRON ultrasonic handpiece.

The I-SURGE micro-motor coupling developed by SATELEC means that the IMPLANT CENTER can be used with most contra-angles (without internal spray or light) found on the market.

The tiltable screen enables the practitioner to work under optimum ergonomic conditions.

In order to get the most out of this high-technology device and ensure it has a long working life, please read this manual carefully before using or maintaining it.

Points preceded by the symbol  should be paid particular attention.

II - WARNINGS

CAUTION:

United States Federal Law restricts the use of this device solely to qualified, trained and competent dental health practitioners or under their supervision.

To reduce the risks of accidents, the precautions stipulated below must be taken.

Device users:

- Use of the IMPLANT CENTER is restricted solely to qualified, trained and competent dental health practitioners in the normal context of their work.
- If you have received this device by error, please contact the supplier so that it can be removed.

Interactions:

- The IMPLANT CENTER must not be used if the patient and/or the operator has a cardiac stimulator (pacemaker) or any other active implant (e.g. a cochlear implant).
- The device complies with applicable electromagnetic compatibility standards. The user should nevertheless ensure that any potential electromagnetic interference does not cause an additional risk (presence of radiofrequency emitters, electronic devices, etc.).
- The device is not designed to withstand shocks delivered by an electric defibrillator.
- Do not attempt to connect to the IMPLANT CENTER any other accessories than those supplied by SATELEC.

Electrical connection:

- Your IMPLANT CENTER must be connected to the electric power supply by a certified dental installation technician.
- To avoid any risk of electric shock, this device must be connected to an electric supply with a protective ground.
- The electric supply to which the device is connected must comply with the standards in force in your country.

Using the device:

- Do not use the device if it appears to be damaged or faulty.
- Turn the device off before unplugging the power cord.
- To unplug the power cord, grip the cord plug and hold the wall socket.
- Never use any other irrigation solution containers other than those intended for suspension from the supplied brackets.
- The device must only be used with bottles or bags of physiological saline.
- The capacity of the irrigation solution containers used must not exceed one litre.
- When the device is not to be used for a long period of time, unplug the device from the electric supply.
- The IMPLANT CENTER must always be used with the screen raised.
- Do not exert excessive force on the tiltable screen.
- Do not move IMPLANT CENTER during use.

Environment:

- Do not cover the device or obstruct the ventilation vents.
- Do not immerse the device in liquid, and do not use it outdoors.
- Do not tilt the device at an angle greater than 5°.
- Do not place the device near a heat source.

- Make sure that the cords are not in a traffic path.
- The device must be stored in its original packaging in an appropriate and safe place.
- The device is not designed for use in the presence of anaesthetic gases or any other flammable gas.
- Do not expose the device to water vapor, or splashes
- The device is not designed to work near ionizing radiation.
- Do not insert metal objects into the device (risk of electric shock, short-circuit or emission of hazardous substances).

Maintenance:

- Before and after each use, your IMPLANT CENTER must be disinfected with products recommended by SATELEC.
- Before each procedure, it is essential to make sure that the accessories to be used have been cleaned, disinfected and sterilized.

Accessories:

- The device can only accept accessories designed by SATELEC for the particular use for which they are intended.
- Use of accessories from other manufacturers is a potential hazard for you or your patients.

Repair:

- Do not repair or modify the device without prior authorization from SATELEC.
- In the case of a fault, contact the supplier of your device. Do not use unauthorized repairers, who might make your device dangerous for you and your patients.

If in any doubt, contact an approved dealer or our customer support department:

- www.acteongroup.com
- E-mail: satelec@acteongroup.com.

III - DESCRIPTION

3. 1 PHYSICAL DESCRIPTION

The IMPLANT CENTER has the following components:

- a control unit (Fig. 1),
- a multifunction footswitch (Fig. 8),
- a micro-motor cord with connectors (Fig. 1-2),
- two brackets for irrigation solution (Fig. 1-3),
- an I-SURGE micro-motor (Fig. 1-4),
- a grounded power cord (Fig. 1-5),
- an LCD (Liquid Crystal Display) screen (Fig. 1-6),
- a control keypad (Fig. 1-7),
- a motor cord connector (Fig. 1-8),
- two peristaltic pump housings (Fig. 1-9),
- a cord connector for the two ultrasonic functions (Fig. 1-10),
- a cord for NEWTRON scaler or PIEZOTOME handpiece (Fig. 1-11) (according to option),
- a NEWTRON handpiece or a PIEZOTOME handpiece with its cord (Fig. 1-12) (according to option).

Located at the back of the unit are:

- 1 power cord socket with ground pin (Figure 2, (1)),
- 1 footswitch connector (Figure 2, (2)),
- 1 fan (Figure 2, (3)),
- 2 bracket holders (Figure 2, (4)),
- 1 power switch (Figure 2, (5)).

3. 2 TECHNICAL DESCRIPTION

a) LCD screen/control keypad

The tiltable control keypad (Fig. 3) is used to define the settings of your IMPLANT CENTER.

The IMPLANT CENTER is adjusted by applying reasonable pressure to the keypad buttons.

Button identification (Fig. 3-1 to 10).

- Item 1: Flush/prime
- Item 2: Increase irrigation flow rate
- Item 3: Reduce irrigation flow rate
- Item 4: Irrigation ON/OFF
- Item 5: Program selection
- Item 6: Save or confirm (depending on mode)
- Item 7: Reduce value
- Item 8: Selection/displacement
- Item 9: Increase value
- Item 10: Change side.

The LCD screen displays the various items of information needed by the user.

Display and indicator identification

Displays common to all modes:

Display and indicator identification (Fig. 4-11 to 15).

- Item 11: Save or confirm (depending on mode)

- Item 12: Irrigation ON/OFF
- Item 13: Error message zone
- Item 14: Irrigation flow rate
- Item 15: Parameter setting zone
- Fig. 3, Item 25: Selected side.

I-SURGE mode:

Display and indicator identification (Fig. 5-16 to 19).

- Item 16: Rotation direction
- Item 17: Control footswitch type
- Item 18: I-SURGE tab
- Item 19: Current program number.

PIEZOTOME mode:

Display and indicator identification (Fig. 6,-20 to 21).

- Item 20: Current program number
- Item 21: PIEZOTOME tab.

NEWTRON mode:

Display and indicator identification (Fig. 7-22 to 234).

- Item 22: Current program number
- Item 23: NEWTRON tab
- Item 24: Color code (CCS) used (Figure 3, (24))

b) Control unit back panel

The power receptacle (Fig. 2-1) with its grounding pin is used to connect the IMPLANT CENTER to the power supply using a plug-in power cord.

The footswitch connector (Fig. 2-2) is used to connect the IMPLANT CENTER to the multifunction control footswitch.

The fan (Fig. 2-3), protected by a metal grille, keeps the IMPLANT CENTER at its optimal performance level.

The bracket holders (Fig. 2-4) are used to install the brackets.

The power switch (Fig. 2-5) is used to switch the device on or off.

c) Control unit sides

- The pump housings (Fig. 1-9) are designed to accommodate SATELEC disposable irrigation line cassettes.
- The irrigation cassettes are installed by lifting the pump covers and inserting the cassettes horizontally into the cassette housings.

d) Front panel

- The left connector is designed to accept the connector of the I-SURGE micro-motor cord. Only a SATELEC micro-motor cord must be used.
- The righthand connector is designed to accept the connector of the NEWTRON cord or the connector of the PIEZOTOME handpiece cord.

e) Control footswitch

The availability of a large number of functions on the control footswitch allows the practitioner to work in a perfectly sterile environment, avoiding the risk of cross-contamination.

Once the various settings have been adjusted, the user has no further need to touch the keypad.

Depending on the mode, the control footswitch buttons have different functions.

I-SURGE mode:

Footswitch button definition (Fig. 8-1 to 5).

- Item 1: Irrigation ON/OFF

- Item 2: Select side
- Item 3: Reverse direction
- Item 4: Change program (1 to 4)
- Item 5: I-SURGE motor control (ON/OFF or progressive)

PIEZOTOME mode:

Footswitch button definition (Fig. 8-1 to 5)

- Item 1: Irrigation ON/OFF
- Item 2: Select side
- Item 3: Flush/prime
- Item 4: Change program (1 to 4)
- Item 5: Ultrasonic ON

NEWTRON mode:

Footswitch button definition (Fig. 8-1 to 5).

- Item 1: Irrigation ON/OFF
- Item 2: Select side
- Item 3: Flush/prime
- Item 4: Change program (1 to 4)
- Item 5: Ultrasonic ON

f) Technical characteristics

Device name: IMPLANT CENTER

Electrical power supply:

Voltage: 100VAC to 230VAC

Frequency: 50Hz/60Hz

Rated power: 250VA @ 230VAC

I-SURGE function

Operation

Intermittent service: 20 sec. ON/30 sec. OFF at 2N.cm

Output characteristics

I-SURGE micro-motor speed: 100rpm to 40,000rpm

I-SURGE micro-motor torque: max. 6N.cm

Irrigation flow rate: 10 to 120ml/min. (nominal value). Adjustment in 5 ml/min. steps

PIEZOTOME function

Operation

Intermittent service: 10min ON/5min OFF

Output characteristics

Voltage: 250 Volts

Min. ultrasonic frequency: 28 kHz

Modulation frequency: 60 Hz/30 Hz/10 Hz according to mode

Irrigation flow rate: 10 to 80 ml/min. (nominal value). Adjustment in 5 ml/min. steps.

NEWTRON function

Operation

Intermittent service: 10 min. ON/5 min. OFF

Output characteristics

No-load voltage: 150 Volts (nominal value without handpiece)

Min. ultrasound frequency: 28 kHz

Irrigation flow rate: 10 to 40 ml/min. (nominal value). Adjustment in 1 ml/min. steps

Protection:

Electrical class: Class 1

Electrical safety class:

Type BF in I-SURGE mode

Type BF in PIEZOTOME mode

Type BF in NEWTRON mode

Safety :

Thermal safety shut off to avoid overheating of the I-SURGE micro-motor.

Safety shut off if internal malfunction.

2 fuses (power cord receptacle):

5 mm x 20 mm/2.5 AT rating 35 A @ 230VAC

1 internal fuse not accessible to the user

Overall dimensions:

Control unit:

Width: 495 mm

Height: 147 mm (LCD screen folded down)

Height with bracket: 444 mm

Depth: 376 mm

Weight: max. 6,000g

Footswitch:

Length: 177 mm

Width: 173 mm

Height: 43 mm

Weight: max. 700 g

LCD screen:

Height: 86 mm

Width: 115 mm

Micro-motor cord: 2,000 mm

Scaler handpiece cord: 2,000 mm

I-SURGE micro-motor:

Length: 99 mm

Diameter: 21 mm

Weight: 120 g

Coupling type: E-type according to standard ISO 3964

Temperatures:

Operating: +10°C to +40°C

Storage: -20°C to +70°C

Humidity:

Operating: 30% to 75%

Storage: 10% to 100% including condensation

Measuring units displayed and meaning:

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rpm = speed (r/min)

g) Materials

The Implant Center keypad and display covering is made of AUTOTEX AM which is a hard-coated, textured polyester film.

IV - INSTALLATION / FIRST USE

4. 1 UNPACKING THE DEVICE

Upon receipt of the device, check for any damage caused in transit.

Contact your supplier if necessary.

4. 2 RECOMMENDATIONS

Have your IMPLANT CENTER connected to the power electricity supply by a certified dental installation technician.

The electrical connection of the IMPLANT CENTER must comply with the applicable standards in your country.

A grounded power supply socket must be used.

4. 3 INSTALLATION

Important:

Do not place the IMPLANT CENTER close to or on top of another device.

Do not place the power cord and the footswitch cord in a wire cover or in a cable gland.

- Place the control unit in position on a fixed surface, horizontal or with a slope of no more than 5°.
- Check that the power switch (Fig. 2-5) is in position 0 (OFF).
- Connect the power cord to the receptacle of the device.
- Connect the power cord to a socket with a ground pin.
- Connect the control footswitch to the footswitch connector (Fig. 2-2).
- Position the footswitch so that it is easily accessible for your feet.
- Insert the brackets into the bracket holders (Fig. 2-4).
- Connect the micro-motor cord to the connector (Fig. 1-8).
- Align the I-SURGE micro-motor and cord connector electrical contacts, then screw the motor onto the connector.
- Connect the NEWTRON scaler cord or the PIEZOTOME handpiece cord to the connector (Fig. 1-10)

(according to option).

- Connect the NEWTRON handpiece to the scaler cord (delivered according to option).
- Place the I-SURGE micro-motor on its rest.
- Place the NEWTRON handpiece on the handpiece support.
- Suspend the bottles or bags of physiological saline from the brackets.
- Raise the LCD screen.
- Adjust the angle of the LCD screen to suit your angle of view.

4. 4 FIRST USE

Before using the IMPLANT CENTER for the first time, all the equipment must be given routine maintenance and/or sterilized according to the procedures defined in chapter 11.

V - IRRIGATION SOLUTION

The IMPLANT CENTER is not designed to administer drugs. It must be used only with bottles or bags of physiological saline or sterile water with an individual capacity not exceeding one liter.

VI - IRRIGATION LINES

The IMPLANT CENTER must be used exclusively with SATELEC irrigation lines.

VII - ADJUSTMENTS / SETTINGS / MODES / INTERFACE

7. 1 STARTUP

By default, the device runs an automatic test sequence (self-check).

Each time the device is switched on, the screen displays program 1 of the I-SURGE function.

7. 2 SETTINGS

The settings are saved when each program is configured and are retrieved whenever that program is selected.

7. 3 IRRIGATION FUNCTION

a) Flush/prime

The flush/prime function is available on the control keypad and on the footswitch.

Pressing button  (Fig. 3-1) or the footswitch button (depending on the function used) switches on the flush/prime function.

The flush/prime function stays on as long as the button is pressed.

b) Flow rate adjustment

The irrigation flow rate is adjusted using the  and  buttons (Fig. 3-2 and 3).

Each time the corresponding buttons are pressed, the flow rate is adjusted.

The adjusted value is displayed on the LCD screen (Fig. 4-14); it is saved in the current program by pressing button .

c) Irrigation ON/OFF

The irrigation is switched on or off by pressing button  or the footswitch button (Fig. 8-1).

When irrigation is selected, the symbol  is displayed.

When irrigation is not selected, the symbol  does not appear on the LCD screen.

7. 4 I-SURGE FUNCTION

a) Selecting the rotation direction

The rotation direction of the I-SURGE micro-motor is selected by the control footswitch (Fig. 8-3).

An audible signal is emitted while the motor is turning in the reverse direction (counterclockwise), and the rotation direction is displayed on the LCD screen: .

b) Program setting adjustment

 **Important:**

When programming the maximum instrument end speed, the entered value cannot exceed the pre-programmed value.

Do not go below the minimum motor speed value (100 rpm, 1:1 contra-angle).

Do not exceed the maximum pre-programmed value of the instrument end torque, and do not go below 10% of this value.

when setting the torque, please refer to the contra-angle and implant system manufacturers' recommendations.

The selected footswitch type applies to all modes. It is not possible to select programs with different footswitch types.

Select the I-SURGE function by pressing the selection button on the left side (Fig. 3-10).

On the LCD screen, the I-SURGE tab must be highlighted (Fig. 5).

The settings are adjusted in the following order:

1- Contra-angle selection

2- Speed adjustment

3- Torque adjustment

4- Footswitch type selection

Select the program for which settings are to be adjusted , ,  or  (Fig. 3-5), then press button ; a cursor appears under CA.

c) Contra-angle selection

Selecting a pre-programmed contra-angle:

Scroll the pre-programmed values using the  or  button (Fig. 3-7 and 9).

The maximum permitted speed and torque for each contra-angle are displayed.

Confirm the selected contra-angle by pressing button .

Programming a user-defined contra-angle:

The IMPLANT CENTER enables specific contra-angles to be used. A contra-angle value can be programmed.

Scroll the preprogrammed values using the  or  button (Fig. 3-7 and 9) until the symbol  appears (Fig. 9-10).

This symbol reminds the user that this is a user-defined mode.

Press button . The cursor is placed under the first figure.

Adjust the value of this figure using the  or  button (Fig. 3-7 and 9).

Press button . The cursor is placed under the second figure.

Adjust the value of the second figure using the  or  button (Fig. 3-7 and 9).

Confirm the final value by pressing button .

Note: While the contra-angle is being adjusted, the maximum permitted speed and torque are displayed automatically.

d) Speed adjustment

When the contra-angle has been selected or entered, the cursor is placed below rpm (Fig. 5). Press button **S**. The cursor is placed under the first figure.

Adjust the value of the first figure using the **+** or **-** button (Fig. 3-7 and 9).

Press button **S**. The cursor is placed under the next figure.

Adjust the speed value figure by using the **+** or **-** and **S** buttons.

Confirm the speed value by pressing button **⇒**.

Note: The instrument speed can be adjusted directly using the **1**, **2** or **3** buttons, even when the multifunction footswitch is switched on.

To obtain a stable speed adjustment result, it is recommended that the footswitch button be fully pressed (Fig. 8-5).

Pressing the **+** or **-** button (Fig. 3-7 and 9) increases or reduces the speed value by 1 rpm.

Pressing and holding the **+** or **-** button (Fig. 3-7 and 9) increases or reduces the speed value.

If the torque or speed value is selected wrongly, !min! or !max! will be displayed, followed up by the maximal or minimal value permitted.

e) Torque adjustment

Press button **S**. The cursor is placed under the first figure.

Adjust the value of the first figure using the **+** or **-** button (Fig. 3-7 and 9).

Press button **S**. The cursor is placed under the next figure.

Repeat this sequence figure by figure.

Confirm the torque value by pressing button **⇒**.

Note: The instrument speed can be adjusted directly using the **+** or **-** button (Fig. 3-7 and 9) for the program **4**, even when the multifunction footswitch is switched on.

f) Footswitch type selection

When the torque has been adjusted, the cursor is displayed under the footswitch symbol.

Press button **S** to toggle between ON/OFF and progressive.

g) Completion of setting adjustment

To be able to use your device with the new settings, press button **⇒**.

The storing of the various settings is confirmed by an audible signal and by the symbol **⇒** which blinks on the screen.

Repeat this sequence if necessary in order to adjust all the programs.

7. 5 PIEZOTOME FUNCTION

⚠ Important:

The **PIEZOTOME** function power settings cannot be adjusted by the user. Only the irrigation settings can be user-defined.

By default, the device is in **PIEZOTOME** mode when no ultrasonic handpiece is connected to it.

The spray for the **PIEZOTOME** function can only be adjusted when a **PIEZOTOME** handpiece is connected to the device.

Select the **PIEZOTOME** function by pressing the right selection button (Fig. 3-10).

On the LCD screen, the **PIEZOTOME** tab must be highlighted (Fig. 6).

Select the program to be used **1**, **2**, **3** or **4** (Fig. 3-5); the selected program is displayed on the LCD screen.

To confirm the choice of program, press button **S**.

If necessary, modify the Irrigation function (flow rate value, etc.) (refer to 6.3).

Confirm the new settings by pressing button **⇒**. The saving of the various settings is confirmed by an audible signal and by the symbol **⇒** which blinks on the screen.

Repeat this sequence if necessary in order to adjust all the programs.

7. 6 NEWTRON FUNCTION

⚠ Important:

A scaler cord and its NEWTRON handpiece must be connected to the device in order to be able to select the various settings.

Select the NEWTRON function by pressing the righthand selection button (Fig. 3-10).

On the LCD screen, the NEWTRON tab must be highlighted (Fig. 7).

Select the program to be adjusted **1**, **2**, **3** or **4** (Fig. 3-5); the selected mode is displayed on the LCD screen and an indicator corresponding to the SATELEC color code (CCS) lights up (Fig. 3-24).

Adjust the power level according to the Tipbook™ recommendations using the **+** or **-** button (Fig. 3-7 and 9). The power level is displayed on the screen (Fig. 7).

If necessary, modify the irrigation function (flow rate value, etc.), then press button **S**.

Confirm the mode settings by pressing button **⇒**. The storing of the various settings is confirmed by an audible signal and by the symbol **⇒** which blinks on the screen.

Repeat this sequence if necessary in order to adjust all the programs.

7. 7 TOOLBOX

⚠ Important:

- The adjusted **TOOLBOX** settings become effective when the device is restarted.
- The **TOOLBOX** function is exited by switching off the device.

The IMPLANT CENTER has a **TOOLBOX** function used to make various adjustments, including contrast, light off time-out (option) and reboot in factory configuration.

To access the **TOOLBOX** function, switch off the device, then press buttons **+** and **-** (Fig. 3-7 and 9) simultaneously while switching the device on again. Release the buttons when the **TOOLBOX** screen is displayed (Fig. 10).

To select a setting, press button **⇒** ; a check mark is displayed next to the selected setting.

a) Contrast adjustment

Confirm access to setting **☞** by pressing button **S**. The cursor is displayed under the contrast value. Adjust the contrast value using the **+** and **-** buttons (Fig. 3-7 and 9).

The contrast is adjustable from 1 (lightest) to 50 (darkest).

Confirm the value by pressing button **⇒**.

b) Light off time-out adjustment

Confirm access to setting **☞** by pressing button **S**. The cursor is displayed under the time-out value. Adjust the time-out value from 0 to 10 by pressing the **+** and **-** buttons (Fig. 3-7 and 9).

The time-out value is adjustable from 0 (no time-out) to 10 (10 seconds).

Confirm the setting by pressing button **⇒**.

c) Reboot in factory configuration

Confirm access to setting  by pressing button  a check mark is displayed next to the selected setting.

Confirm the return to the factory configuration by pressing button .

d) Software version

The version of the embedded software (X.Y) is displayed next to the symbol .

The IMPLANT CENTER is configured in the factory with the following settings:

Motor mode:

Program	Function	Contra-angle ratio	Speed at instrument end	Torque at instrument end	Irrigation
1	Marking	20 :1	1200 tr/min	80 N.cm	80 ml/min
2	Pilot drilling	20 :1	800 tr/min	80 N.cm	100 ml/min
3	Tapping	20 :1	15 tr/min	40 N.cm	100 ml/min
4	Screwing	20 :1	30 tr/min	40 N.cm	0

PIEZOTOME mode:

Program	Function	Modulation frequency	Irrigation
1	D1	60 Hz	50 ml/min
2	D2	30 Hz	50 ml/min
3	D3	10 Hz	50 ml/min
4	D4	30 Hz	50 ml/min

NEWTRON mode:

Program	Function		Power	Irrigation
1	Green	Soft	P = 5	15 ml/min
2	Yellow	Medium	P = 5	15 ml/min
3	Blue	High	P = 5	15 ml/min
4	Orange	Boost	P = 5	15 ml/min

VIII - SAFETY

The IMPLANT CENTER has a system that detects malfunctions of the device.

When an excessively high temperature of the I-SURGE micro-motor is detected, symbol 1 in Fig. 9 is displayed in the error zone (Fig. 4-13) and the device beeps three times.

The device switches to degraded mode so that it can complete the dental procedure.

The available torque value is set at 25% in order to protect the I-SURGE micro-motor.

It is recommended that the I-SURGE micro-motor be allowed to cool down until symbol 1 in Fig. 9 disappears.

When an internal operation error occurs, symbol 2 in Fig. 9 is displayed in the error zone and the device beeps three times.

It is recommended that the device be switched off (Fig. 2-5) and switched back on again.

IX - USING THE DEVICE

Important:

- Do not disconnect the micro-motor cord or the scaler cord when the device is switched on and the footswitch is pressed.
- Do not disconnect the micro-motor or the ultrasonic handpieces when the device is switched on and the footswitch is pressed.
- Do not engage the instrument to or disengage it from the contra-angle when the micro-motor is running.
- Do not screw or unscrew tips when the handpieces are on.
- All accessories must be cleaned, disinfected and sterilized before use.
- For your safety and that of your patient, the IMPLANT CENTER must not be used with accessories other than those supplied or recommended by SATELEC.
- Check the integrity of the device and its accessories before and after each use in order to detect any problem.
- If necessary, do not use the device and replace any faulty component.

The IMPLANT CENTER is prepared for use as follows:

- Make sure that the device is connected correctly and properly insulated.
- Switch the power switch (Fig. 2-5) to position I (power on).
- Suspend bottles or bags of physiological saline from the brackets.
- Raise the irrigation pump covers.
- Insert the irrigation line cassettes horizontally into their housings.
- Close the covers.
- Insert the irrigation line perforating pins into the bottles or pouches of physiological saline.
- Raise the LCD screen.

- Adjust the angle of the LCD screen according to the position of your device and your angle of view.
- Make any necessary setting adjustments (irrigation flow rate, speed, torque, power, power mode, program, etc.) (refer to chapter 7).

Note: This sequence must be completed under conditions of controlled asepsis for the safety of your patients. An operating assistant may be needed.

Note: the data is saved automatically, in order to maintain the display parameters as previously used upon switchover from the various modes.

I-SURGE function

- Check that there are no traces of moisture on the I-SURGE micro-motor or micro-motor cord connections. Remove any moisture by wiping and by blowing using the multifunction syringe fed with medical quality filtered air.
- Connect the micro-motor cord to the left connector of the device (Fig. 1-8).
- Align the micro-motor and cord connector electrical contacts, then screw the motor onto the connector.
- Connect a contra-angle to the I-SURGE micro-motor.
- Fit an instrument (drill, bur, etc.) to the contra-angle.
- Fasten the irrigation line to the cord using line clips.
- Connect the end of the irrigation line to the contra-angle.
- Run the motor fitted with its contra-angle and its instrument without irrigation at moderate speed for 10 to 15 seconds in order to spread the lubricant and eliminate excess.
- Prime the irrigation circuit by pressing the control keypad flush button (Fig. 3-1).
- When the irrigation solution reaches the end of

the contra-angle, release the footswitch.

- Use the IMPLANT CENTER in accordance with good dental practice.

PIEZOTOME function

- Connect the PIEZOTOME handpiece cord to the right connector of the device (Fig. 1-10).
- Screw the selected tip onto the handpiece using a wrench (according to the recommendations of the clinical booklets).
- Fasten the irrigation line to the cord using line clips.
- Connect the end of the irrigation line to the PIEZOTOME handpiece.
- Prime the irrigation circuit by pressing the flush button on the multifunction footswitch.
- When the irrigation solution reaches the end of the tip, release the footswitch.
- Use the IMPLANT CENTER in accordance with good dental practice.

NEWTRON function

- Check that there are no traces of moisture on the NEWTRON handpiece or the scaler cord connections. Remove any moisture by wiping and by blowing using the multifunction syringe fed with medical quality filtered air.
- Connect the scaler cord to the right connector of the device (Fig. 1-10).
- Align the electrical contacts of the NEWTRON handpiece and the scaler cord, then plug the handpiece into the cord.
- Screw the selected tip onto the handpiece using a wrench (according to the Tipbook recommendations).
- Fasten the irrigation line to the cord using line clips.
- Connect the end of the irrigation line to the scaler cord.
- Prime the irrigation circuit by pressing the Flush button on the multifunction footswitch.
- When the irrigation solution reaches the end of

the tip, release the footswitch.

- Use the IMPLANT CENTER in accordance with good dental practice.

X - SHUTTING DOWN THE DEVICE

Upon completion of the dental procedure:

- Remove the bottles or bags of physiological saline from the brackets.
- Remove the irrigation line perforating pins from the bottles or bags.
- Immerse the irrigation line perforating pins in a vessel containing distilled water.
- Rinse each irrigation line alternately, along with the contra-angle and the handpiece, by running the flush function until the vessel and the irrigation lines are completely empty.
- Remove the irrigation line clips.
- Disconnect the single-use irrigation lines from the contra-angle and from the handpiece and dispose of them in a safety container for soiled medical items.
- Remove the rotary instruments fastened to the contra-angle and the tip screwed onto the handpiece.
- Disconnect the contra-angle from the I-SURGE micro-motor.
- Disconnect the micro-motor cord from the I-SURGE.
- Disconnect the PIEZOTOME handpiece or disconnect the NEWTRON handpiece and the scaler cord.
- Switch the device off (0).

XI - ROUTINE MAINTENANCE / STERILIZATION

Important:

Except for the irrigation lines, no accessory is delivered sterile.

During sterilization, items made of different metals must not come into contact. Any contact would establish galvanic couples and cause localized damage.

To avoid this happening, place the items in individual sterilizable bags or in a sterilization box.

To maintain the sterility or asepsis of the accessories (contra-angle, micro-motor, micro-motor cord, handpieces, etc.), store them in sealed bags or containers suitable for medical practice.

The maintenance and/or sterilization instructions below must be applied before each use of the device.

Avoid using cleaning and disinfection agents containing flammable substances.

If this is not possible, make sure that all such agent has evaporated and that there is no combustible matter on the device or its accessories before it is switched on.

Sterilized items must be allowed to dry and cool to room temperature before they are reused.

Before sterilization, check that your autoclave is clean and that the water used is of adequate quality.

After each sterilization cycle, remove the items from the autoclave immediately in order to minimize the risks of metal corrosion.

11. 1 IRRIGATION LINE MAINTENANCE

The irrigation lines supplied by SATELEC are disposable and must always be discarded after use.

Re-use of an irrigation line may result in contamination of your patients.

Do not attempt to resterilize the single-use irrigation lines.

Do not attempt to modify the irrigation lines.

11. 2 CONTRA-ANGLE MAINTENANCE

Refer to the manufacturer's instructions for your contra-angle.

11. 3 ROTARY INSTRUMENT MAINTENANCE

Please refer to the manufacturer's instructions for your instruments (burs, drills, etc.).

11. 4 DEVICE MAINTENANCE

Important:

Never use an abrasive cleaning agent on the device.

Do not allow cleaning/disinfecting spray to penetrate inside the device.

The device, the bracket and the multifunction footswitch are not sterilizable.

The IMPLANT CENTER control unit, the control footswitch and the bracket must always be cleaned and disinfected after each procedure, using alcohol, disinfectant agents or disinfectant wipes for dental practices, such as SEPTOL™ wipes.

Regular monitoring of the IMPLANT CENTER unit is necessary in order to detect any problem. It is important to keep the control unit ventilation vents clean in order to avoid abnormal heating.

11. 5 MICRO-MOTOR CORD AND SCALER CORD MAINTENANCE

Important:

Never use an abrasive cleaning agent on the cords.

The cords must not be immersed.

The cords must not be disinfected using agents containing acetone, chlorine or bleach.

The cords must not be cleaned in an ultrasonic bath.

a) Cleaning and disinfection

The cords must always be cleaned, disinfected and sterilized after each procedure, using alcohol, disinfectant agents or disinfectant wipes for dental practices, such as SEPTOL™ wipes.

b) Sterilization:

The cords can be sterilized in an autoclave under the following conditions:

- Autoclave: Class B.
- Sterilization temperature: 134°C at 2 bars.
- Time at sterilization temperature: 18 minutes.

11. 6 I-SURGE MICRO-MOTOR MAINTENANCE

Important:

Never use an abrasive cleaning agent on the I-SURGE micro-motor.

The I-SURGE micro-motor must not be immersed or disinfected using agents containing acetone, chlorine or bleach.

The I-SURGE micro-motor must not be cleaned in an ultrasonic bath or disassembled.

The I-SURGE micro-motor is a high-technology device. To maintain optimal performance, the following recommendations must be complied with.

a) Cleaning and disinfection

Use cleaning and disinfection agents compatible with the I-SURGE micro-motor (such as SPRAYNET®). When using a spray, spray the outside and inside of the I-SURGE micro-motor moderately to expel residues.

Clean and disinfect the surfaces using a clean, lint-free cloth impregnated with the agent.

b) Lubrication

Lubricate before each sterilization or at least twice a day.

Use lubricants compatible with the I-SURGE micro-motor (such as LUBRIFLUID®).

Place a clean cloth over the I-SURGE micro-motor ventilation vents to absorb the surplus lubricant. Inject the lubricant through the nose of the micro-motor using an oilcan of lubricant (or a spray) for one-half to one second.

c) Sterilization

The I-SURGE micro-motor can be sterilized in an autoclave under the following conditions:

- Autoclave: Class B.
- Sterilization temperature: 134°C.
- Time at sterilization temperature: 18 minutes.

11. 7 ULTRASONIC HANDPIECE MAINTENANCE

After each use, the handpiece irrigation circuit must be rinsed with distilled or demineralized water for 20 to 30 seconds.

Cleaning and disinfection:

The handpieces must always be cleaned, disinfected and sterilized after each procedure, using alcohol, disinfectant agents or disinfectant wipes for dental practices, such as Septol wipes.

For handpiece sterilization, refer to the manual for each product.

11. 8 TIP MAINTENANCE

a) Lifetime

Tip shape and weight are the determining characteristics for obtaining maximal efficiency of the ultrasonic generator.

Careful monitoring of these two characteristics by the user will ensure that optimal performance of the device is maintained.

Consequently, you are strongly advised not to alter the structure of the tips by filing or twisting them or modifying them in any other way.

Similarly, through normal wear, tip aging leads to alteration of its characteristics.

Always replace a tip that has been damaged by wear or by accidental impact (dropped, bent, etc.).

Use the Tipcard™ to check scaler tip wear. Replace routinely-used tips at least once to three times a year.

b) Pre-disinfection/chemical cleaning

To be carried out immediately after treatment, wearing thick gloves.

Immerse in an ultrasonic bath (enzymatic, quaternary ammonium base detergent/disinfectant solution), carefully respecting the concentration and contact time recommended by the manufacturer of the solution.

Use a solution that bears the EC marking or that is in compliance with any standard required by the national regulations.

Rinse under running water for at least 30 seconds.

c) Mechanical/chemical cleaning

To be carried out immediately after pre-disinfection/chemical cleaning, wearing thick gloves.

Scrub the products in a new bath of enzymatic, quaternary ammonium base detergent/disinfectant solution, for at least 30 seconds, and in any case until all remaining traces of contamination have been removed.

Rinse again under running water for at least 30 seconds.

d) Drying

Dry with a clean, non-woven single-use cloth, so as to remove any traces of liquid.

Pack in single-use sterilization pouches or bags in compliance with the specifications defined in standard EN ISO 11607-1 or in any equivalent standard required by the national regulations.

e) Sterilization

The tips must be sterilized individually in an autoclave, according to the following parameters:

- autoclave: type B, compliant with standard EN 13060;
- sterilization temperature: 134°C;

- sterilization stabilization time: 18 minutes;
- pressure: 2 Bars min.

f) Storage

Then, store the sterilized products in a dry place protected from dust.

Before re-use, if the integrity of the packaging is not certain, package then re-sterilize in accordance with the defined protocol.

If contamination is visible inside the bag, destroy the product.

g) Disposal of the product

Dispose of the product in receptacles for waste materials of healthcare activities involving infectious risks.

CAUTION: Not all autoclaves can reach 134°C. Not all autoclaves draw a pre-vacuum. Please refer to your autoclave manufacturer for specific sterilization instructions.

XII - MONITORING/ PREVENTIVE AND CORRECTIVE MAINTENANCE

Important:

In the case of a fault, you are advised to contact the supplier of your device. Do not use unauthorized repairers, who might make your device dangerous for you and your patients.

12. 1 MONITORING

Regular monitoring of the device and its accessories is necessary in order to detect any insulation fault or damage. Replace them if necessary.

It is important to keep the control unit ventilation vents clean in order to avoid abnormal heating.

12. 2 PREVENTIVE AND CORRECTIVE MAINTENANCE

SATELEC recommends inspection or overhaul of the I-SURGE micro-motor at least once a year.

12. 3 FUSE REPLACEMENT

Your IMPLANT CENTER is protected by two fuses located in the power cord (Fig. 2-1).

The replacement procedure is as follows:

- Switch the device off (position O).
- Unplug the power cord from the electricity supply.
- Unplug the power cord from the receptacle (Fig. 2-1).
- Insert the tip of a flat screwdriver into the notch above the fuse drawer to release it.
- Remove the old fuses.
- Replace the old fuses with fuses of the same type and rating.
- Push the fuse drawer back into its housing until it clicks to indicate the correct position.
- Plug the power cord into the receptacle (Fig. 2-1).
- Plug the power cord into the electrical supply.

Note: The device also has an internal fuse not accessible to the user.

Contact SATELEC after-sales service (see chapter 2 - Repairs).

All the information necessary for the repair of faulty components is made available by SATELEC to the technical personnel of the SATELEC-approved dealer network who carry out the work.

12. 4 OPERATION FAULTS

Refer to the tables on the next page.

Fault detected	Possible causes	Solutions
Device does not operate (LCD screen off)	Power cord not plugged in properly.	- Check the electric socket. - Return to SATELEC after-sales.
	Power switch in position 0.	Place the power switch in position I.
	No power voltage.	Call an electrician.
	Power receptacle fuse(s) blown.	Replace the power receptacle fuses.
	Internal fuse blown.	Return to SATELEC® after-sales.
Device does not operate (LCD screen on)	Transmission fault.	- Switch the device off then on again. - Return to SATELEC after-sales.
Torque fault	Thermal cut-out.	Allow the micro-motor to cool down.
	Incorrect torque adjustment.	Adjust the torque in accordance with good dental practice.
	Wrong contra-angle.	- Change the contra-angle. - Return to SATELEC® after-sales.
Speed fault	Incorrect speed adjustment.	Adjust the speed in accordance with good dental practice.
	Wrong contra-angle.	- Change the contra-angle. - Return to SATELEC after-sales.
No spray	Irrigation solution bag or bottle empty.	Replace the bag or bottle of irrigation solution.
	Irrigation off.	Press the irrigation ON/OFF button.
	Irrigation line clogged.	Change irrigation line.
	Line end piece on contra-angle clogged.	Clear the contra-angle line end piece.
Unsuitable spray	Incorrect irrigation flow rate adjustment.	Adjust the irrigation flow rate.

Fault detected	Possible causes	Solutions
Micro-motor does not operate	Faulty micro-motor cord connector electrical contacts.	Clean the electrical contacts of the micro-motor cord connectors.
	Faulty micro-motor electrical contacts.	Clean the electrical contacts of the micro-motor.
	Micro-motor cord wires cut.	Return to SATELEC after-sales to replace the cord.
	Motor damaged.	Return to SATELEC after-sales to replace the motor.
Power fault weak tip vibration.	Tip worn or distorted.	Replace the tip.
	Incorrect power adjustment.	Refer to the Tipbook for the necessary information.
	Incorrect use: incorrect angle of attack or inadequate pressure on the tooth.	Refer to the Tipbook for the necessary information.
	Liquid or moisture present between handpiece and cord.	Dry the electrical contacts thoroughly.
No ultrasonic function	Incorrect tip tightening.	Retighten the tip with the wrench.
	Faulty connector contact.	Clean the connector contacts.
	Handpiece cord wire cut.	Return to SATELEC after-sales to replace the cord.
Leak between the NEWTRON handpiece and the handpiece cord	Worn handpiece 1.15x1 seal.	Replace the seal (kit F12304).
Liquid leak in the irrigation pump	Broken tube in the irrigation line cassette.	Replace the irrigation line.

XIII - ELECTROMAGNETIC COMPATIBILITY

Important:

The power cord, the micro-motor cord, the ultrasonic handpiece cords and the multifunction footswitch cord must be kept away from each other.

The IMPLANT CENTER requires special precautions to be taken with regard to electromagnetic compatibility.

It must be installed and prepared for use as described in chapter 4.

Certain types of mobile telecommunication devices such as mobile telephones are likely to interfere with the IMPLANT CENTER.

The recommended separation distances in this paragraph must therefore be complied with.
The IMPLANT CENTER must not be used near or on top of another device.
If this cannot be avoided, its operation under the conditions of use must be checked beforehand.
The use of accessories other than those specified or sold by Satelec® as replacement parts may have the consequence of increasing the emissions or decreasing the immunity of the IMPLANT CENTER.

13. 1 ELECTROMAGNETIC EMISSIONS

The IMPLANT CENTER is intended for use in the electromagnetic environment specified in the table below.
The user and/or installer must ensure that the IMPLANT CENTER is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF emission - CISPR 11.	Group 1	IMPLANT CENTER uses RF energy for internal operation. Therefore, its radiofrequency emissions are very low and are not likely to cause any interference in nearby equipment.
RF emission - CISPR 11.	Class B	IMPLANT CENTER is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic current emission IEC 61000-3-2.	Class A	
Voltage fluctuation and flicker IEC 61000-3-3.	Complies	

13. 2 ELECTROMAGNETIC IMMUNITY

The IMPLANT CENTER is intended for use in the electromagnetic environment specified in the table below.
The user and/or installer must ensure that the IMPLANT CENTER is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2.	± 6KV contact ± 8KV air	± 6KV contact ± 8KV air	Floors must be wood, concrete, cement or tiled. If floors are covered with synthetic material (carpet, etc.), the relative humidity must be at least 30%.
Electrical fast transients IEC 61000-4-4.	± 2KV for power supply lines	± 2KV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment (hospital, clinic).
Surges IEC 61000-4-5.	± 1KV differential mode ± 2KV common mode	± 1KV differential mode ± 2KV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations IEC 61000-4-11.	<5% UT (>95% dip in UT) for 0.5 cycles. 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 250 cycles	<5% UT (>95% dip in UT) for 0.5 cycles. 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 250 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the use of the IMPLANT CENTER requires continued operation during power mains interruption, it is recommended that the product be powered from a separate power supply (UPS, etc.).

13. 3 ELECTROMAGNETIC IMMUNITY / RADIOFREQUENCY MOBILE EQUIPMENT

The IMPLANT CENTER is intended for use in the electromagnetic environment specified in the table below. The user and/or installer must ensure that the IMPLANT CENTER is used in such an electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Portable and mobile radiofrequency communications devices must not be used near the IMPLANT CENTER (including its cables) at a distance less than that recommended and calculated according to the frequency and power of the emitter.			
Conducted disturbance, radiofrequency fields. IEC61000-4-6	3 V/m 150 KHz to 80 MHz	3 V/m	Recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated radiofrequency electromagnetic field. IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80MHz à 800MHz. $d = 2.3 \sqrt{P}$ 800MHz à 2.5GHz. Where P is the maximum power rating of the emitter in watts (W) according to the manufacturer's specifications and d is the recommended minimum separation distance in meters (m).
The electromagnetic field strengths of fixed radiofrequency emitters, as determined by an electromagnetic environment measurement (a), must be less than the compliance level in each frequency range (b). Interference may occur near equipment marked with the symbol below: 			

- Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.
- Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

- (a): The electromagnetic field strengths of fixed radiofrequency emitters, such as base stations for mobile telephones (cellular/cordless), mobile radios, amateur radio, AM/FM radio broadcasts and TV broadcasts cannot be determined exactly by theory. To assess the electromagnetic environment due to fixed radiofrequency emitters, an electromagnetic environment measurement must be made. If the measured radiofrequency field strength in the immediate environment where the product is used exceeds the compliance level specified above, the performance of the product must be tested to verify whether it conforms to the specification. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the product.
- (b): In the 150kHz to 80MHz frequency range, the electromagnetic field strengths must be less than 3V/m.

13. 4 RECOMMENDED SEPARATION DISTANCES

The IMPLANT CENTER is intended for use in an electromagnetic environment in which radiated radiofrequency disturbances are controlled.

The IMPLANT CENTER user and/or installer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile radiofrequency communications equipment emitters and the IMPLANT CENTER, according to the maximum output power of the equipment, as recommended in the table below.

Rated max. power of the emitter (W)	Separation distance in metres (m) according to emitter frequency		
	150 KHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0,01	0.12 m	0.12 m	0.23 m
0,1	0.38 m	0.38 m	0.73 m
1	1.2 m	1.2 m	2.3 m
10	3.8 m	3.8 m	7.3 m
100	12 m	12 m	23 m

For emitters rated at max. power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the emitter, where P is the max. power rating of the emitter in watts (W) according to the manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

10. 5 CABLE LENGTHS

Cables and accessories	Maximum length	Complies with:
Handpiece cords	Less than 3m	RF emission, CISPR 1 - Class B/Group 1
Footswitch cord		Harmonic current emission: IEC61000-3-2
Power supply cord		Voltage fluctuation: IEC61000-3-3
		Immunity to electrostatic discharge: IEC61000-4-2
		Immunity to electrical fast transients/bursts: IEC61000-4-4
		Immunity to surges: IEC61000-4-5
		Immunity to voltage dips, short interruptions and voltage variations: IEC61000-4-11
	Immunity to conducted disturbances induced by radiofrequency fields: IEC61000-4-6	
	Immunity to radiated radiofrequency electromagnetic fields: IEC61000-4-3	

XIV - DISPOSAL AND RECYCLING

As electrical and electronic equipment, the device must be disposed of according to a specialized procedure for collection, pick-up and recycling or destruction (in particular on the European market, with reference to Directive 2002/96/EC of 23/01/2003).

When your device reaches the end of its life, we consequently recommend that you contact your dental equipment dealer (or, failing this, the nearest ACTEON Group office, the list of which is given in chapter 18), for information on how to proceed.

XV - MANUFACTURER'S LIABILITY

This equipment has been designed and developed according to an ISO 13485-certified quality assurance system.

The manufacturer is not liable if:

- the manufacturer's installation recommendations have not been followed (supply voltage, electromagnetic environment, etc.),
- repairs or other work have been done on the device by persons not authorized by the manufacturer,
- the device has been used connected to an electrical system that does not comply with current regulations,
- the device has been used in ways other than those specified in this manual,
- accessories (tips, handpiece, irrigation lines, etc.) other than those supplied by SATELEC have been used,
- the instructions in this document have not been followed.

Note: The manufacturer reserves the right to modify the device and/or the user's manual without notice.

XVI - REGULATIONS

This medical device is classified as class IIa according to European Directive 93/42/EEC.

This equipment is manufactured in compliance with the current IEC 60601-1 standard.

XVII - SYMBOLS

	Alternating current
	BF device
	Warning, please refer to the accompanying documentation
	Not to be disposed of as household waste
	CE marking

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